

hIL10(19-178)-His

Cat. No.

Product information


Source:	HEK293
Extinction Coefficient:	0.395
AA MW:	19kD
Concentration:	lot specific
Purity:	>95% determined by reducing SDS-PAGE >90% determined by SEC-HPLC
Endotoxin Level:	<10EU/mg determined by LAL method
Sterility:	0.2 μ M filtered

Formulation and Storage


Formulation:	PBS, pH 7.4
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: DAD1A, Sig=280, 4 Ref=off

RT[min]	Width[min]	Area	Height	Area%
6.433	1.47	481.76	19.64	23.13
7.033	1.46	1600.69	53.93	76.87
SUM	2.93	2082.45		

- For Research Use Only. Not for use in diagnostic procedures.