

Protein name His-TEV-Kras(2-185) G12D

Cat. No.

Product information

Nature Recombinant

Expression system: Escherichia coli

Concentration: lot specific

Purity: >95% determined by SDS-PAGE
>95% determined by SEC-HPLC

Species Human

Extinction Coefficient: 13660

AA MW: 23319.34Da

Preparation and Storage


Format: Aqueous buffer solution

Formulation: 20 mM HEPES pH 7.4, 150 mM NaCl, 0.5 mM MgCl₂, 1 mM TCEP


Storage and handling: Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: DADIA, Sig=280, 4 Ref=off

RT [min]	Width [min]	Area	Height	Area%
6.897	2.85	11113.66	383.87	100.00
SUM	2.85	11113.66		

- For Research Use Only. Not for use in diagnostic procedures.