

**Protein name** His-TEV-Kras(2-185) G12V

**Cat. No.**

### Product information


**Nature** Recombinant  
**Expression system:** Escherichia coli  
**Concentration:** lot specific  
**Purity:** >90% determined by SDS-PAGE  
**Species** Human  
**Extinction Coefficient:** 13660  
**AA MW:** 23303.38Da

### Preparation and Storage


**Format:** Aqueous buffer solution  
**Formulation:** 20 mM HEPES pH 7.4, 150 mM NaCl, 0.5 mM MgCl<sub>2</sub>, 1 mM TCEP  
**Storage and handling:** Store at -80°C. Avoid multiple freeze-thaw cycles.

### Assay Data

## SDS-PAGE



## SEC-HPLC



Signal: DAD1A, Sig=280, 4 Ref=off

RT[min]	Width[min]	Area	Height	Area%
6.897	4.39	15751.32	547.68	100.00
SUM	4.39	15751.32		

- For Research Use Only. Not for use in diagnostic procedures.