

Protein name Flag-TEV-SOS1(564-1094)

Cat. No.

Product information


Nature	Recombinant
Expression system	Escherichia coli
Concentration:	lot specific
Purity:	>95% determined by SDS-PAGE >95% determined by SEC-HPLC
Extinction Coefficient:	54570
AA MW:	61990.94Da

Formulation and Storage


Format:	Aqueous buffer solution
Formulation:	20 mM Tris, pH 8.0, 100 mM NaCl, 10%Glycerol
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: DAD1A, Sig=280, 4 Ref=off

RT[min]	Width[min]	Area	Height	Area%
5.601	2.58	5266.32	179.97	100.00
SUM	2.58	5266.32		

- For Research Use Only. Not for use in diagnostic procedures.