

Benralizumab

Cat. No.

Product information


Isotype:	hIgG1
Source:	CHO-K1
Extinction Coefficient:	1.44
AA MW:	146kD
Concentration:	lot specific
Purity:	>95% determined by SDS-PAGE >95% determined by SEC-HPLC
Endotoxin Level:	<1EU/mg determined by LAL method
Sterility:	0.2 μ M filtered

Formulation and Storage


Formulation:	PBS, pH 7.4
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: VWD1A, Wavelength=280 nm

RT[min]	Width[min]	Area	Height	Area%
4.197	1.21	75.27	2.51	1.32
5.202	3.30	5628.85	188.08	98.68
SUM	4.51	5704.12		

- For Research Use Only. Not for use in diagnostic procedures.

Bioactivity