

Dupilumab

Cat. No.

Product information


Isotype:	hIgG4
Source:	CHO-K1
Extinction Coefficient:	1.412
AA MW:	147kD
Concentration:	lot specific
Purity:	>95% determined by SDS-PAGE >95% determined by SEC-HPLC
Endotoxin Level:	<1EU/mg determined by LAL method
Sterility:	0.2 μ M filtered

Formulation and Storage


Formulation:	PBS, pH 7.4
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: WD1A, Wavelength=280 nm

RT[min]	Width(min)	Area	Height	Area%
4.324	1.09	165.97	5.80	3.92
5.029	3.35	4068.69	149.35	96.08
SUM	4.44	4234.66		

- For Research Use Only. Not for use in diagnostic procedures.

Bioactivity