

Mirikizumab

Cat. No.

Product information


Isotype:	hIgG4
Source:	CHO-K1
Extinction Coefficient:	1.533
AA MW:	144kD
Concentration:	lot specific
Purity:	>95% determined by SDS-PAGE >95% determined by SEC-HPLC
Endotoxin Level:	<1EU/mg determined by LAL method
Sterility:	0.2 μ M filtered

Formulation and Storage


Formulation:	PBS, pH 7.4
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: VWD1A, Wavelength=280 nm

RT[min]	Width[min]	Area	Height	Area%
4.306	0.94	24.31	1.38	0.56
4.965	3.54	4300.78	122.60	99.44
SUM	4.48	4325.09		

- For Research Use Only. Not for use in diagnostic procedures.

Bioactivity