

Omalizumab

Cat. No.

Product information


Isotype:	hIgG1
Source:	CHO-K1
Extinction Coefficient:	1.565
AA MW:	147kD
Concentration:	lot specific
Purity:	>95% determined by SDS-PAGE >95% determined by SEC-HPLC
Endotoxin Level:	<1EU/mg determined by LAL method
Sterility:	0.2µM filtered

Formulation and Storage


Formulation:	PBS, pH 7.4
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: YND1A, Wavelength=280 nm

RT[min]	Width[min]	Area	Height	Area%
4.164	0.54	11.94	0.56	0.19
4.873	3.50	6224.75	245.02	99.81
SUM	4.04	6236.70		

- For Research Use Only. Not for use in diagnostic procedures.

Bioactivity