

Plozalizumab

Cat. No.

Product information


Isotype:	hIgG1
Source:	CHO-K1
Extinction Coefficient:	1.414
AA MW:	146kD
Concentration:	lot specific
Purity:	>95% determined by SDS-PAGE >95% determined by SEC-HPLC
Endotoxin Level:	<1EU/mg determined by LAL method
Sterility:	0.2µM filtered

Formulation and Storage


Formulation:	PBS, pH 7.4
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: VRDIA, Wavelength=280 nm

RT[min]	Width[min]	Area	Height	Area%
4.469	0.58	19.43	1.17	0.30
4.917	2.60	6515.81	291.90	99.70
SUM	3.18	6535.24		

- For Research Use Only. Not for use in diagnostic procedures.

Bioactivity