

Ranibizumab

Cat. No.

Product information


Isotype:	hKappa
Source:	CHO-K1
Extinction Coefficient:	1.806
AA MW:	48kD
Concentration:	lot specific
Purity:	>95% determined by SDS-PAGE >95% determined by SEC-HPLC
Endotoxin Level:	<1EU/mg determined by LAL method
Sterility:	0.2 μ M filtered

Formulation and Storage


Formulation:	PBS, pH 7.4
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: YVDIA, Wavelength=280 nm

RT[min]	Width[mm]	Area	Height	Area%
6.524	0.66	176.92	10.26	2.58
7.011	2.07	6692.77	214.39	97.42
SUM	2.73	6869.69		

- For Research Use Only. Not for use in diagnostic procedures.

Bioactivity