

hCD19(20-291)-His

Cat. No.

Product information


Source:	HEK293
Extinction Coefficient:	2.516
AA MW:	31kD
Concentration:	lot specific
Purity:	>90% determined by reducing SDS-PAGE polymers determined by SEC-HPLC
Endotoxin Level:	<10EU/mg determined by LAL method
Sterility:	0.2µM filtered

Formulation and Storage


Formulation:	PBS, pH 7.4
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: DMDA, Sig=280, 4 Ref=off

RT[min]	Width[min]	Area	Height	Area%
4.730	1.70	5905.28	109.08	55.34
5.032	2.61	4765.85	113.42	44.66
SUM	4.31	10671.13		

- For Research Use Only. Not for use in diagnostic procedures.