

Vedolizumab

Cat. No.

Product information


Isotype:	hIgG1
Source:	CHO-K1
Extinction Coefficient:	1.562
AA MW:	147kD
Concentration:	lot specific
Purity:	>95% determined by SDS-PAGE >95% determined by SEC-HPLC
Endotoxin Level:	<1EU/mg determined by LAL method
Sterility:	0.2 μ M filtered

Formulation and Storage


Formulation:	PBS, pH 7.4
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: DAD1A, Sig=280,4 Ref=off

Retention Time[min]	Peak Width[min]	Peak Area	Peak Height	Peak Area%	Name
4.392	0.74	43.86	1.86	0.78	
5.168	2.10	5602.78	212.43	99.22	
	SUM	5646.64			

- For Research Use Only. Not for use in diagnostic procedures.

Bioactivity