

Nirsevimab

Cat. No.

Product information


Isotype:	hIgG1
Source:	CHO-K1
Extinction Coefficient:	1.376
AA MW:	146kD
Concentration:	lot specific
Purity:	>95% determined by SDS-PAGE >95% determined by SEC-HPLC
Endotoxin Level:	<1EU/mg determined by LAL method
Sterility:	0.2 μ M filtered

Formulation and Storage


Formulation:	PBS, pH 7.4
Storage and handling:	Store at -80°C. Avoid multiple freeze-thaw cycles.

Assay Data

SDS-PAGE



SEC-HPLC



Signal: VMD1A, Wavelength=280 nm

RT[min]	Width[min]	Area	Height	Area%
4.242	0.76	39.99	1.76	0.89
4.781	3.52	4431.98	134.61	99.11
SUM	4.28	4471.97		

- For Research Use Only. Not for use in diagnostic procedures.

Bioactivity